US/CA/UK/FR/ES/PT/DE/IT/RO

IN260300230V01

713-221V00

04 STEPS

10MIN

EN_ **IMPORTANT, RETAIN FOR FUTURE REFERENCE: READ CAREFULLY.**

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EN

WARNING

DO NOT use if poorly attached. If you are unable to securely fasten all parts, contact a dealer or trained technician to assemble. Take care during assembly not to lose your balance when testing the connections.
DO NOT use if any parts were damaged in shipping or become worn from use. Contact your dealer for replacement and repair before further use.
ALWAYS use the toilet seat when sitting during use. Sitting off-center or directly on the bucket may cause damage or personal injury.
DO NOT use on uneven ground or with legs at uneven heights. All four legs' rubber tips MUST touch the floor at all times.
ALWAYS take care to remain balanced while sitting and rising, using both armrests.
DO NOT exceed this product's weight limit.
NOT for use by children. Keep away from children to reduce the risk of damage or injury.

MAINTENANCE AND CARE

Clean & dry the bucket promptly after each use.
After shower, please wipe the chair clean and keep it dry.
Store in a dry-well ventilated place away from direct sunlight.

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④ x01

Done

EN_ Place the main frame upright, then insert backrest ④ into the main frame.

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US

If you have any questions, please contact our customer care center, Our contact details are below:

001-877-644-9366

customerservice@aosom.com

Imported by Aosom LLC
7777 Center Ave. Suite 430
Huntington Beach, CA 92647-9998 USA
MADE IN CHINA

CA

If you have any questions, please contact our customer care center, Our contact details are below:

416-792-6088

customerservice@aosom.ca

Imported by Aosom Canada Inc.
7270 Woodbine Avenue, Unit 307, Markham,
Ontario Canada L3R 4B9
MADE IN CHINA

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UK

If you have any questions, please contact our customer care center, Our contact details are below:

0044-800-240-4004

enquiries@mhstar.co.uk

IMPORTER ADDRESS:
MH STAR UK Ltd
1 Northampton Cross Logistics Park
Northampton Cross Road
Northampton
NN4 9FH
MADE IN CHINA

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MH Handel GmbH